
2CI30/2CI31/2CI31-P

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Contenido :

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

2CI30/2CI31/2CI31-P Assembling Instructions

Contents :

- 2 blue joints
- 2 yellow dummies
- 4 nuts
- 4 screws
- 4 nails
- 1 washer to support the bushing
- 1 bumper

Place the first dummy on the lateral malleolus and fix it with nails. Position the second one facing it on the opposite side of the positive model and fix it with nails.

Fill with plaster the empty areas that may be left between the dummy and the positive model.

To create the housing for the extension bumper, place on the positive model a clean and warm block (170-190 °C) of the same material as the one used for thermoforming, so that the total thickness and the total width exceed 10 mm. The block height must allow a 20-25 mm extra-thickness after cutting off.

We recommend polypropylene PL140xx or PL111xx for thermoforming (xx = material thickness).

In case of thermoforming with Plastazote (PE foam), raise the dummy. Carry out vacuum thermoforming.

After cooling, trace the orthosis cutting lines and cut off plastic with the dummies remaining at their place.


Strip the boot and remove the dummy residues. Polish the edges according to the cutting lines.

Take the lower part of the orthosis and drill (Ø 3.1) the extra-thickness on 10 mm. Bevel with a Ø 8 drill. Stick the bumper.

Drill the clearance holes (Ø 4) for the screws where the dummy has left depressions. Place the nuts and position the joint.

Position the screws and tighten them using an hex wrench and the washer used to support the bushing.

Secure the screws with Loctite ref. XC047 during the final assembly.

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