hyTREK

**GENOU HYDRAULIQUE  
MONOAXIAL  
1P130 / 1P130-KD**

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
hyTREK

**SINGLE AXIS  
HYDRAULIC KNEE  
1P130 / 1P130-KD**

*Instructions for Use by the Prosthetist*

*Protective hood to be glued if necessary*

**150 kg / 330 lbs**

*Adjustment wrench (included)*

## 1 – Description and functions

1P13099-0318 6/20

This knee is intended exclusively for lower limb prosthesis.

It is a single axis hydraulic knee with hydraulic braking in stance phase, which provides 4 different adjustments, an extension assist and a locking system for flexion.

The 4 adjustments are indicated by the following letters :

- R** - adjusts resistance to flexion intensity during stance phase
- S** - adjusts the threshold value for entering into swing phase (sensitivity)
- F** - adjusts flexion angle during swing phase
- E** - adjusts extension dampening at terminal impact

## 2 – Indications, maximum weight

It is specifically recommended for patients with medium to very high activity level (K3/K4) up to 150 Kg / 330 lbs (carried load included).

This knee weighs 1265 g for 1P130 model and 1285 g for 1P130-KD model.

Maximum knee flexion : 120°. However it can be limited by the socket dimensions or by the cosmetic cover.

## 3 – Parts for connection and finishing (Fig. 1a and 1b)

Possible connections :

This knee is available in two different models :

- **1P130** : Upper connection by a male pyramid, convenient for most transfemoral amputation or hip disarticulation.
- **1P130-KD** : Threaded upper connection (M36×1.5), convenient for long residual limb or knee disarticulation.

With these two models, a 34mm Ø tube is used for the lower connection (See Fig. 1b for tightening - See Fig.2 for restrictions).

RELATED PRODUCTS RECOMMENDED :

	1P130	1P130-KD
<b>Upper connection</b>		
Lamination adapters	1K179+1K183 (max weight 150Kg) 1K173 (maximum weight 125Kg) 1K176 (maximum weight 125Kg) 1K177 (maximum weight 100Kg)	1K179 (maximum weight 150Kg)
Adapters	1K172 (maximum weight 150Kg) 1K207-HD (maximum weight 150Kg)	1K184 (maximum weight 150Kg)
<b>Lower connection</b>		
Tube and adapter	1G01-HD (maximum weight 150Kg) 1D52-P6 (maximum weight 125Kg)	
Feet	Dynatrek - 1A600 (maximum weight 150Kg) Dynastar - 1A500 (maximum weight 125Kg) Dyna C - 1A400 (maximum weight 125Kg) DynaCity - 1A510 (maximum weight 125Kg)	

**FINISHING:** One-part cosmetic cover 1G16 (complete or partial)

It is recommended to reinforce the cosmetic cover in front of the joint to increase its product life.

If the prosthesis is worn without cosmetic cover, the supplied protective hood can be glued with XC050.

The cosmetic cover can modify the initial adjustments of the knee and limit some motions. In any case you should place the cosmetic cover before checking the knee functions.

Do not remove the protective devices of the knee, especially the rear band that protects the hydraulic system during full flexion.

## 4 – Alignments (Fig.2)

**Static alignments :**

In a sagittal plane, the load line goes through the Great Trochanter, from 0 to 5 mm ahead of the knee axis. At the foot level, comply with the manufacturer instructions for alignment.

In a frontal plane, the load line goes through the socket center, the knee center and between the first and the second toes.

**Dynamic alignements :**


Bringing the knee forwards decreases stability of the assembly and eases step completion during walking.

Bringing the knee backwards increases stability and the patient's feeling of safety during walking, but makes it harder to start flexion.

The position of the adapter on the socket is important. If needed, use a translation device during fitting trial in order to find the best alignment.

The knee adjustments may require an adjustment of the dynamic alignments (see § 5).


**Fig. 1a**


**Fig. 1b**

**Order of screw tightening :**

- 1) Tighten screw V1 up to 5 Nm
- 2) Tighten screw V2 up to 5 Nm
- 3) Retighten screw V1 up to 5 Nm

**Note :** The M36x1.5 threaded connection of the KD model does not allow the alignment to be modified.  
Check the alignment with an appropriate tool (laser, plumb-line, ...)

## 5 – Adjustments (Fig.3)

This knee allows 4 possible adjustments.

To ensure the patient's safety, the prosthetist should perform the first fitting trials between parallel bars and with factory settings (See Fig.3). Then make the adjustments preferably in the following order:

1. Adjustment R of resistance to flexion intensity during stance phase
2. Adjustment S of the threshold value for entering into swing phase (sensitivity)
3. Adjustment F of flexion angle during swing phase
4. Adjustment E of extension dampening at terminal impact

The foot model used can have an impact on the knee adjustments. The knee should be adjusted each time the foot is changed. The use of energy return feet such as Dynatrek, Dyna C, Dynastar or DynaCity is recommended.

Ensure the patient's safety in all circumstances and especially during slope or stair descent. A new adjustment of the knee is necessary after each maintenance operation.

An 2,5 mm hex wrench comes with the knee and allows to make all the adjustments.

### 1. Adjustment of resistance intensity during stance phase (screw R)

The 1P130 knee comes with a fairly high adjustment of resistance intensity (Fig. 3). The first adjustments should be made carefully and very safely to avoid any risk of falling for the patient.

Screw in by successive 1/12 turn to increase resistance during stance phase (braking). Unscrew to do the opposite.

**Adjustment range = 3/4 turn (please note : 1 turn = same adjustment)**

Between parallel bars, have the patient sit down on a chair in order to set a first braking resistance. In a second time after the other adjustments, refine the resistance adjustment during stair descent with alternate steps, then during slope descent while ensuring the patient's safety with a handrail.

## 2. Adjustment of threshold value for entering into swing phase (screw S)

1P13099-0318 8/20

 The 1P130 knee comes with factory settings that are convenient for most of the users (Fig. 3). The dynamic alignment should allow transition into swing phase and braking activation. The threshold can be modified if the swing phase does not come at toes off, or if resistance does not appear at heel strike.

Screw in by successive 1/4 turn until the swing phase gets released at toes off. Unscrew by successive 1/4 turn to ease brake activation.

**Adjustment range = -1 turn / +2 turns**

 Do not fully loosen screw **S** because this would damage the mechanism. If in doubt, reset the adjustment so that the screw head be as flush as possible.  
There is a real risk of falling if this adjustment is wrong.

## 3. Adjustment of flexion angle during swing phase (screw F) :

Screw in by successive 1/12 turn to limit flexion angle during swing phase, especially with high walking speed. Unscrew to do the opposite.

**Adjustment range = 3/4 turn**


(please note : 1 turn = same adjustment)

## 4. Adjustment of extension dampening at terminal impact (screw E) :

Screw in by successive 1/4 turn to increase dampening at terminal impact (end of extension). Unscrew to do the opposite.

**Adjustment range = ±1 turn**

Once these adjustments are performed, make sure that the knee can return to full extension at low walking speed to avoid any risk of falling for the patient.

## Maximal flexion limit (Fig.4)

In full flexion, the socket can touch the posterior band located above the adjustment knobs.

 The socket should not touch the hydraulic system to preserve it and to avoid knee damages.

The knee is equipped with an inside stop that limits flexion to 120°.

## Lock (Fig.4)

The 1P130 knee is equipped with a lock that can be activated when the user is standing or sitting.

To lock the knee, press knob  at the back of the knee. Always check carefully this operation before using the knee. Press the opposite knob to release the knee.

## 6 – Product life

This component has been tested in accordance with standard ISO10328, load level P7 (150Kg / 330 lbs) ; the cyclical test of 3 million cycles corresponds to a use of 4 or 5 years depending on the patient activity.

## 7 – Maintenance

A yearly control of the knee is recommended to check proper functioning and to clean it (dust).

If necessary, proceed to a new adjustment of the knee.

Flexion stops (RL0402) and protective hood (1P13055) can be replaced using the appropriate cyanoacrylate glue (XC050).

## 8 – Advice for use, maintenance and safety

- In order not to damage the knee, do not use talcum powder to eliminate possible rubbing noise, but rather a silicon spray. Talcum powder damages the mechanical parts, which may lead to function defects and risk of falling for the patient.

**The use of talcum powder will render all claims against PROTEOR null and void.**


- It is strictly forbidden to remove the protective covers and to tighten or loosen any screw of this knee except for the 4 adjustment screws **R, S, F, and E** (Fig.3) and the 2 tightening screws **V1 and V2** (Fig. 1b) of the Ø34 tube holder.
- Never grease the knee axis and ball bearings because this could damage them rapidly.

### **Inform the patient :**


- That the warranty is not applicable to damages due to improper use, inadequate alignment, use in a very dusty environment and without adapted protection, or any inappropriate use.
- That there is a risk of having a finger caught or clothing pinched in the knee. To avoid any risk of injury caused by the moving joint, take care that no one puts his finger near or inside the mechanism.
- That it should be avoided to expose the knee in environments that can cause corrosion of metal parts (fresh water, sea water, chlorinated water, acids, etc.). The knee is resistant to rain but needs to be dried after exposition to water.
- When the patient is carrying a load, the knee functions can be disturbed. Especially resistance during stair descent can be insufficient or the knee can get stuck suddenly.
- That in case of fall, shock or abnormal behaviour of the prosthesis, the knee should be checked by the prosthetist.
- That holding on to the handrail is recommended when going downstairs to avoid risks of falling.
- That at first use after a rest period of several hours, the knee can emit a clicking noise.
- That the temperature range for knee use is from -10 °C to + 40 °C (from 14°F to 104°F).

## 9 – Recycling - End of life

This product includes items made of various materials : elastomer, plastics, aluminum, titanium, steel and brass. It also contains oil. They must be recycled in compliance with the laws in force.

**Fig. 4 :**  
Back view  
of the knee

**Fig. 5**

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